

SPRADLING VALENCIA VINYL SELECTION

Flammability Rating: BS5852 Ignition Source (Crib) 5

Anti-Microbial Protection: Anti-fungal, mildew resistant ASTM-G21. Anti-bacterial, AATCC Test Test Method 147. Permablok 3 from Spradling is a proven vinyl protective coating that's engineered to create a tough, effective barrier against the three biggest problems encountered in hospital and healthcare environments - germs, abrasion and stains.

Abrasion Resistance: >300,000 Martindale Rubs.

Cleaning: Wipe clean routinely with a microfibre cloth, lightly dampened with cold water only, then dry with a clean, absorbent cloth. Remove body fluids and extreme soiling immediately (within 15 minutes) with cold water. Apply your approved infection protocol, a 5% sodium chlorite bleach solution, Haz-Tabs or Chlor-Clean. Remove product residues and then dry with a clean, absorbent cloth.

Stains: Dyes and pigments from jeans, ball point and felt tip pens may permanently stain.

Colours: Although every effort is made to ensure that vinyls supplied match the colours on this page, variations may occur due to limitations in printing processes.

